

Product Name: TP8L1 Rabbit Polyclonal Antibody
Catalog #: APRab19145



Summary

Production Name	TP8L1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	TNFAIP8L1
Alternative Names	
Gene ID	126282.0
SwissProt ID	Q8WVP5. Synthesized peptide derived from human protein . at AA range: 40-120

Application

Dilution Ratio	WB 1:500-2000, ELISA 1:5000-20000
Molecular Weight	20kDa

Background


similarity:Belongs to the TNFAIP8 family.,similarity:Belongs to the TNFAIP8 family.,

Product Name: TP8L1 Rabbit Polyclonal Antibody
Catalog #: APRab19145



Research Area

Image Data



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

Note

For research use only.