

Product Name: KD-Validated UBE2F Recombinant Rabbit Monoclonal Antibody**Catalog #: KVA00366**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:200-1:2,000; ICC 1:100-1:1,000
Molecular Weight	Calculated MW: 21.1kDa

Antigen Information

Gene Name	UBE2F UBE2F; Ubiquitin Conjugating Enzyme E2 F (Putative); NCE2; Ubiquitin-Conjugating Enzyme E2F (Putative); RING-Type E3 NEDD8 Transferase UBE2F; NEDD8-Conjugating Enzyme
Alternative Names	UBE2F; NEDD8 Carrier Protein UBE2F; NEDD8 Protein Ligase UBE2F; NEDD8-Conjugating Enzyme 2; Ubiquitin-Conjugating Enzyme E2 F; NEDD8 Conjugating Enzyme; EC 2.3.2.34; EC 6.3.2
Gene ID	140739.0
SwissProt ID	Q969M7
Immunogen	Recombinant protein of human UBE2F

Background

Enables NEDD8 conjugating enzyme activity. Involved in protein neddylation. Predicted to be active in cytosol and nucleus.
[provided by Alliance of Genome Resources, Jul 2025]

Research Area

Epigenetics and Nuclear Signaling, Cell Biology

Image Data

