

Product Name: KD-Validated HNRNPUL1 Recombinant Rabbit Monoclonal Antibody
Catalog #: KVA00487

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:200-1:2,000
Molecular Weight	Calculated MW: 95.7kDa

Antigen Information

Gene Name	HNRNPUL1
Alternative Names	HNRNPUL1; Heterogeneous Nuclear Ribonucleoprotein U Like 1; E1B-AP5; E1BAP5;
	HNRPUL1; Heterogeneous Nuclear Ribonucleoprotein U-Like Protein 1; Adenovirus Early Region 1B-Associated Protein 5; E1B-55 KDa-Associated Protein 5; E1B 55kDa Associated Protein 5; FLJ12944
Gene ID	11100.0
SwissProt ID	Q9BUJ2
Immunogen	A synthesized peptide derived from human E1B-AP5

Background

This gene encodes a nuclear RNA-binding protein of the heterogeneous nuclear ribonucleoprotein (hnRNP) family. This protein

binds specifically to adenovirus early-1B-55kDa oncoprotein. It may play an important role in nucleocytoplasmic RNA transport, and its function is modulated by early-1B-55kDa in adenovirus-infected cells. [provided by RefSeq, Mar 2016]

Research Area

Cell Biology

Image Data

