

Product Name: KD-Validated HNRNPA1L2 Recombinant Rabbit Monoclonal Antibody
Catalog #: KVAb00495

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:100-1:1,000; ICC 1:100-1:1,000
Molecular Weight	Calculated MW: 34.2kDa

Antigen Information

Gene Name	HNRNPA1L2
Alternative Names	HNRNPA1L2; Heterogeneous Nuclear Ribonucleoprotein A1 Like 2; Heterogeneous Nuclear Ribonucleoprotein A1-Like 2; HnRNP Core Protein A1-Like 2; HnRNP A1-Like 2; LOC144983; HNRNPA1L
Gene ID	144983.0
SwissProt ID	Q32P51
Immunogen	A synthesized peptide derived from human Heterogeneous Nuclear Ribonucleoprotein (A1-like)

Background





Predicted to enable RNA binding activity. Predicted to be involved in mRNA splicing, via spliceosome. Predicted to be located


in cytoplasm and nucleus. Predicted to be part of catalytic step 2 spliceosome. [provided by Alliance of Genome Resources, Jul 2025]

Research Area

Signal Transduction

Image Data





Flow cytometric analysis of HNRNPA1L2 expression in HepG2 cells using HNRNPA1L2 antibody (KVA00495, 1:1,000). Green, isotype control; red, HNRNPA1L2.