
Product Name: KD-Validated CCDC93 Mouse Monoclonal Antibody**Catalog #: KVA00584**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Mouse
Application	WB,FCM
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Mouse mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:100-1:1,000
Molecular Weight	Calculated MW: 73.2kDa

Antigen Information

Gene Name	CCDC93
Alternative Names	CCDC93; Coiled-Coil Domain Containing 93; Coiled-Coil Domain-Containing Protein 93; FLJ10996
Gene ID	54520.0
SwissProt ID	Q567U6
Immunogen	Recombinant protein of human CCDC93


Background

Involved in Golgi to plasma membrane transport and endocytic recycling. Located in intracellular membrane-bounded organelle. [provided by Alliance of Genome Resources, Apr 2025]


Research Area

Cell Biology


Image Data



Western blotting analysis using coiled-coil domain containing 93 antibody (KVA00584). Coiled-coil domain containing 93 expression in wild type (WT) and coiled-coil domain containing 93 (CCDC93) shRNA knockdown (KD) HT-1080 cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with coiled-coil domain containing 93 antibody (KVA00584, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (APS0631, 1:10,000) respectively.



Western blotting analysis using coiled-coil domain containing 93 antibody (KVA00584). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with coiled-coil domain containing 93 antibody (KVA00584, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (APS0631, 1:10,000) respectively.



Flow cytometric analysis of Coiled-coil domain containing 93 expression in HeLa cells using Coiled-coil domain containing 93 antibody (KVA00584, 1:1,000). Green, isotype control; red, Coiled-coil domain containing 93.

Coiled-coil domain containing 93-Alexa Fluor® 647