

Product Name: KD-Validated ERBB2 Mouse Monoclonal Antibody**Catalog #: KVA00692**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Mouse
Application	WB,FCM
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Mouse mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:2,500-1:5,000; FC 1:200-1:2,000
Molecular Weight	Calculated MW: 137.9kDa

Antigen Information

Gene Name	ERBB2
Alternative Names	ERBB2; Erb-B2 Receptor Tyrosine Kinase 2; HER2; NEU 2; C-ERB-2; C-ERB2; MLN-19; HER-2; CD340; NGL; V-Erb-B2 Avian Erythroblastic Leukemia Viral Oncogene Homolog 2; Tyrosine Kinase-Type Cell Surface Receptor HER2; Neuro/Glioblastoma Derived Oncogene Homolog; Human Epidermal Growth Factor Receptor 2; Receptor Tyrosine-Protein Kinase ErbB-2; Metastatic Lymph Node Gene 19 Protein; Proto-Oncogene C-ErbB-2; Proto-Oncogene Neu; P185(ErbB2); EC 2.7.10.1; MLN 19; V-Erb-B2 Avian Erythroblastic Leukemia Viral Oncogene Homolog 2 (Neuro/Glioblastoma Derived Oncogene Homolog) 2; V-Erb-B2 Erythroblastic Leukemia Viral Oncogene Homolog 2, Neuro/Glioblastoma Derived Oncogene Homolog; V-Erb-B2 Avian Erythroblastic Leukemia Viral Oncoprotein 2; Neuroblastoma/Glioblastoma Derived Oncogene Homolog; Metastatic Lymph Node Gene 19; C-Erb B2/Neu Protein;

CD340 Antigen; P185(ERBB2); HER-2/Neu; Herstatin; P185erbB2; EC 2.7.10; VSCN2; MLN19; TKR1

Gene ID	2064.0
SwissProt ID	P04626
Immunogen	Recombinant protein of human ERBB2


Background

This gene encodes a member of the epidermal growth factor (EGF) receptor family of receptor tyrosine kinases. This protein has no ligand binding domain of its own and therefore cannot bind growth factors. However, it does bind tightly to other ligand-bound EGF receptor family members to form a heterodimer, stabilizing ligand binding and enhancing kinase-mediated activation of downstream signalling pathways, such as those involving mitogen-activated protein kinase and phosphatidylinositol-3 kinase. Allelic variations at amino acid positions 654 and 655 of isoform a (positions 624 and 625 of isoform b) have been reported, with the most common allele, Ile654/Ile655, shown here. Amplification and/or overexpression of this gene has been reported in numerous cancers, including breast and ovarian tumors. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized. [provided by RefSeq, Jul 2008]


Research Area

Tags & Cell Markers, Signal Transduction, Cancer, Neuroscience


Image Data



Western blotting analysis using erb-b2 receptor tyrosine kinase 2 antibody (KVA00692). Erb-b2 receptor tyrosine kinase 2 expression in wild type (WT) and erb-b2 receptor tyrosine kinase 2 (ERBB2) shRNA knockdown (KD) 293T cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with erb-b2 receptor tyrosine kinase 2 antibody (KVA00692, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (APS0631, 1:10,000) respectively.




Western blotting analysis using erb-b2 receptor tyrosine kinase 2 antibody (KVA00692). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with erb-b2 receptor tyrosine kinase 2 antibody (KVA00692, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (APS0631, 1:10,000) respectively.



Validation of Erb-b2 receptor tyrosine kinase 2 knockdown using flow cytometry. Wild-type (WT, Blue) and knockdown (KD, Green) 293T cells were stained with Erb-b2 receptor tyrosine kinase 2 antibody (KVA00692, 1:2,000) and analyzed using BD flow cytometer.

Erb-b2 receptor tyrosine kinase 2-Alexa Fluor® 647



Flow cytometric analysis of Erb-b2 receptor tyrosine kinase 2 expression in H9c2 cells using Erb-b2 receptor tyrosine kinase 2 antibody (KVA00692, 1:2,000). Green, isotype control; red, Erb-b2 receptor tyrosine kinase 2.