

Product Name: KD-Validated TATA Element Modulatory Factor 1 Recombinant Rabbit Monoclonal Antibody

Catalog #: KVAAb00839

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:200-1:2,000; ICC 1:100-1:1,000
Molecular Weight	Calculated MW: 122.8kDa

Antigen Information

Gene Name	TMF1 TMF1; TATA Element Modulatory Factor 1; ARA160; TMF; Androgen Receptor-Associated
Alternative Names	Protein Of 160 KDa; Androgen Receptor Coactivator 160 KDa Protein; TATA Element Modulatory Factor
Gene ID	7110.0
SwissProt ID	P82094
Immunogen	A synthesized peptide derived from human TMF

Background





Enables nuclear androgen receptor binding activity and transcription coactivator activity. Involved in androgen receptor


signaling pathway and positive regulation of transcription by RNA polymerase II. Located in Golgi apparatus. [provided by Alliance of Genome Resources, Jul 2025]

Research Area

Signal Transduction, Epigenetics and Nuclear Signaling,

Image Data





Flow cytometric analysis of TMF expression in HepG2 cells using TMF antibody (KVA00839, 1:2,000). Green, isotype control; red, TMF.