

Product Name: KD-Validated RABEPK Recombinant Rabbit Monoclonal Antibody
Catalog #: KVA00911

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:200-1:2,000
Molecular Weight	Calculated MW: 40.6kDa

Antigen Information

Gene Name	RABEPK
Alternative Names	Rab9 Effector Protein With Kelch Motifs; RAB9P40; BA65N13.1; 40 KDa Rab9 Effector Protein; P40; Rab9 Effector P40
Gene ID	10244.0
SwissProt ID	Q7Z6M1
Immunogen	A synthesized peptide derived from human p40

Background

Predicted to be involved in receptor-mediated endocytosis and vesicle docking involved in exocytosis. Located in intracellular membrane-bounded organelle. [provided by Alliance of Genome Resources, Apr 2025]

Research Area

Signal Transduction, Cell Biology

Image Data

