

Product Name: KD-Validated Phospho-Paxillin (Y118) Recombinant Rabbit Monoclonal Antibody

Catalog #: KVA00966

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000
Molecular Weight	Calculated MW: 64.5kDa

Antigen Information

Gene Name	PXN
Alternative Names	Paxillin; Testicular Tissue Protein Li 134
Gene ID	5829.0
SwissProt ID	P49023
Immunogen	A synthesized peptide derived from human Phospho-Paxillin (Y118)

Background


This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties

(PMID:9054445). [provided by RefSeq, Aug 2011]


Research Area

Signal Transduction,Cancer

Image Data



Western blotting analysis using Phospho-Paxillin (Y118) antibody (KVA00966). Phospho-Paxillin (Y118) expression in wild type (WT) and Phospho-Paxillin (Y118) shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. β -Tubulin serves as a loading control. The blot was incubated with Phospho-Paxillin (Y118) antibody (KVA00966, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (APS0635, 1:10,000) respectively.



Western blotting analysis using Phospho-Paxillin (Y118) antibody (KVA00966). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with Phospho-Paxillin (Y118) antibody (KVA00966, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (APS0635, 1:10,000) respectively.