
Product Name: KD-Validated Caspase 6 Recombinant Rabbit Monoclonal Antibody
Catalog #: KVAb01037

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:200-1:2,000; ICC 1:100-1:1,000
Molecular Weight	Calculated MW: 33.3kDa

Antigen Information

Gene Name	CASP6
Alternative Names	CASP6; Caspase 6; CSP-6; MCH2; Caspase-6; Caspase 6, Apoptosis-Related Cysteine Peptidase; Caspase 6, Apoptosis-Related Cysteine Protease; Mammalian Ced-3 Homologue; EC 3.4.22.59; Apoptotic Protease MCH-2; Apoptotic Protease Mch-2; CASPASE-6; EC 3.4.22;
	CASP-6
Gene ID	839.0
SwissProt ID	P55212
Immunogen	A synthesized peptide derived from human Caspase-6

Background


This gene encodes a member of the cysteine-aspartic acid protease (caspase) family of enzymes. Sequential activation of

caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic acid residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein is processed by caspases 7, 8 and 10, and is thought to function as a downstream enzyme in the caspase activation cascade. Alternative splicing of this gene results in multiple transcript variants that encode different isoforms. [provided by RefSeq, Oct 2015]


Research Area

Neuroscience, Cell Biology


Image Data




Western blotting analysis using caspase 6 antibody (KVA01037). Caspase 6 expression in wild-type (WT) and caspase 6 (CASP6) shRNA knockdown (KD) C2C12 cells with 20 μ g of total cell lysates. β -Tubulin serves as a loading control. The blot was incubated with caspase 6 antibody (KVA01037, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (APS0635, 1:10,000) respectively.




Western blotting analysis using Caspase 6 antibody (KVA01037). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with Caspase 6 antibody (KVA01037, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (APS0635, 1:10,000) respectively.



Immunocytochemical staining of HT-1080 cells with Caspase 6 antibody (KVA01037, 1:1,000). Nuclei were stained blue with DAPI; Caspase 6 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser intensity and smart gain: Medium. Scale bar: 20 μ m.



Validation of Caspase 6 knockdown using flow cytometry. Wild-type (WT, Blue) and knockdown (KD, Green) HeLa cells were stained with Caspase 6 antibody (KVA01037, 1:2,000) and analyzed using BD flow cytometer.



Flow cytometric analysis of Caspase 6 expression in HT-1080 cells using Caspase 6 antibody (KVA01037, 1:2,000). Green, isotype control; red, Caspase 6.