
Product Name: KD-Validated PPWD1 Mouse Monoclonal Antibody**Catalog #: KVA01132**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Mouse
Application	WB,FCM
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Mouse mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2,500; FC 1:200-1:2,000
Molecular Weight	Calculated MW: 73.6kDa

Antigen Information

Gene Name	PPWD1 PPWD1; Peptidylprolyl Isomerase Domain And WD Repeat Containing 1; KIAA0073;
Alternative Names	Peptidylprolyl Isomerase Domain And WD Repeat-Containing Protein 1; Spliceosome-Associated Cyclophilin; EC 5.2.1.8; Epididymis Secretory Sperm Binding Protein
Gene ID	23398.0
SwissProt ID	Q96BP3
Immunogen	Recombinant protein of human PPWD1

Background

Enables peptidyl-prolyl cis-trans isomerase activity. Predicted to be involved in protein folding. Located in nuclear body. Part of catalytic step 2 spliceosome. [provided by Alliance of Genome Resources, Jul 2025]

Research Area

Epigenetics and Nuclear Signaling

Image Data

