

Product Name: KD-Validated TRIP13 Mouse Monoclonal Antibody**Catalog #: KVA01228**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Mouse
Application	WB,FCM
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Mouse mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2,500; FC 1:200-1:2,000
Molecular Weight	Calculated MW: 48.6kDa

Antigen Information

Gene Name	TRIP13 TRIP13; Thyroid Hormone Receptor Interactor 13; 16E1BP; Human Papillomavirus Type 16 E1 Protein-Binding Protein; Thyroid Receptor-Interacting Protein 13; Pachytene Checkpoint
Alternative Names	Protein 2 Homolog; TR-Interacting Protein 13; 16E1-BP; TRIP-13; Thyroid Receptor Interacting Protein 13; HPV16 E1 Protein Binding Protein; HPV16 E1 Protein-Binding Protein; SH3TC1/TRIP13 Fusion; OZEMA9; OOMD9; MVA3; PCH2
Gene ID	9319.0
SwissProt ID	Q15645
Immunogen	Recombinant protein of human TRIP13

Background

This gene encodes a protein that interacts with thyroid hormone receptors, also known as hormone-dependent transcription factors. The gene product interacts specifically with the ligand binding domain. This gene is one of several that may play a role in early-stage non-small cell lung cancer. [provided by RefSeq, Oct 2009]

Research Area

Neuroscience

Image Data

