
Product Name: KD-Validated RAB5C Recombinant Rabbit Monoclonal Antibody**Catalog #: KVA01341**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:200-1:2,000
Molecular Weight	Calculated MW: 23.5kDa

Antigen Information

Gene Name	RAB5C RAB5C; RAB5C, Member RAS Oncogene Family; RAB5CL; RABL; RAB5C, Member Of RAS
Alternative Names	Oncogene Family; Ras-Related Protein Rab-5C; L1880; RAB5L; RAB, Member Of RAS Oncogene Family-Like; EC 3.6.5.2
Gene ID	5878.0
SwissProt ID	P51148
Immunogen	A synthesized peptide derived from human RAB5C / RABL

Background

Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment (Han et al., 1996 [PubMed 8646882]).[supplied by

OMIM, Nov 2010]

Research Area

Signal Transduction

Image Data

