

Product Name: KD-Validated CD97 Recombinant Rabbit Monoclonal Antibody**Catalog #: KVA01535**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Rabbit
Application	WB,FCM,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rabbit IgG
Clonality	Rabbit mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1,000-1:5,000; FC 1:200-1:2,000; ICC 1:100-1:1,000
Molecular Weight	Calculated MW: 91.9kDa

Antigen Information

Gene Name	ADGRE5
Alternative Names	ADGRE5; Adhesion G Protein-Coupled Receptor E5; Leukocyte Antigen CD97; TM7LN1; CD97; Seven-Transmembrane, Heterodimeric Receptor Associated With Inflammation; Seven Transmembrane Helix Receptor; Seven-Span Transmembrane Protein; CD97 Molecule; CD97 Antigen
Gene ID	976.0
SwissProt ID	P48960
Immunogen	A synthesized peptide derived from human CD97

Background


This gene encodes a member of the EGF-TM7 subfamily of adhesion G protein-coupled receptors, which mediate cell-cell

interactions. These proteins are cleaved by self-catalytic proteolysis into a large extracellular subunit and seven-span transmembrane subunit, which associate at the cell surface as a receptor complex. The encoded protein may play a role in cell adhesion as well as leukocyte recruitment, activation and migration, and contains multiple extracellular EGF-like repeats which mediate binding to chondroitin sulfate and the cell surface complement regulatory protein CD55. Expression of this gene may play a role in the progression of several types of cancer. Alternatively spliced transcript variants encoding multiple isoforms with 3 to 5 EGF-like repeats have been observed for this gene. This gene is found in a cluster with other EGF-TM7 genes on the short arm of chromosome 19. [provided by RefSeq, Jun 2011]


Research Area

Neuroscience


Image Data




Western blotting analysis using CD97 antibody (KVA01535). CD97 expression in wild type (WT) and CD97 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with CD97 antibody (KVA01535, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (APS0635, 1:10,000) respectively.



Western blotting analysis using CD97 antibody (KVA01535). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with CD97 antibody (KVA01535, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (APS0635, 1:10,000) respectively.



Immunocytochemical staining of HepG2 cells with CD97 antibody (KVA01535, 1:1,000). Nuclei were stained blue with DAPI; CD97 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser intensity and smart gain: Medium. Scale bar: 20 µm.



Flow cytometric analysis of CD97 expression in HepG2 cells using CD97 antibody (KVA01535, 1:2,000). Green, isotype control; red, CD97.