

Product Name: KD-Validated EPS15L1 Mouse Monoclonal Antibody**Catalog #: KVA01942**

For research use only.

Summary

Description	KO&KD-Validated antibody
Host	Mouse
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Mouse mAb
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:400-1:2,000
Molecular Weight	Calculated MW: 94.3kDa

Antigen Information

Gene Name	EPS15L1
Alternative Names	EPS15L1; Epidermal Growth Factor Receptor Pathway Substrate 15 Like 1; Eps15R; Epidermal Growth Factor Receptor Substrate 15-Like 1; Eps15-Related Protein; Epidermal Growth Factor Receptor Pathway Substrate 15-Like 1; Epidermal Growth Factor Receptor Substrate EPS15R; EPS15R
Gene ID	58513.0
SwissProt ID	Q9UBC2
Immunogen	Recombinant protein of human EPS15L1

Background

Enables cadherin binding activity. Predicted to be involved in endosomal transport and synaptic vesicle endocytosis. Located in

membrane. [provided by Alliance of Genome Resources, Apr 2025]

Research Area

Signal Transduction

Image Data

