

Magnetic Separation Rack

Catalog No.: BRC0026

Product Introduction

Magnet is a magnetic separation rack for immune cell sorting, also known as magneticrack, magnetic sorting rack, etc., which can be used with EnkiLife's magnetic bead separation products. It is used for 5mL flow cytometry tubes, Single Well.

Storage Conditions

Storage Conditions: Store at ambient temperature of 15-25°C, Please keep away from pacemakers, magnets, electronic products such as computers and watches, and other objects that respond to magnetic fields.

Shipping: RT

Precautions

1. For magnetic cell separation racks, EnkiLife's magnetic rack (Cat. No.: BRC0026) is the preferred recommendation. STEMCELL's magnetic rack is also acceptable. If the customer does not have access to either, please inform them that the magnetic field strength of the rack must be >7000Gs; otherwise, cell separation efficiency will be compromised.
2. EnkiLife uses magnetic beads with a particle size of 300nm, and therefore is not compatible with Miltenyi Biotec brand magnetic racks.