

Product Name: IL-2 Receptor alpha Mouse Monoclonal Antibody
Catalog #: AMM00727



Summary

Production Name	IL-2 Receptor alpha Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	IL2RA
Alternative Names	IL2RA; Interleukin-2 receptor subunit alpha; IL-2 receptor subunit alpha; IL-2-RA; IL-2R subunit alpha; IL2-RA; TAC antigen; p55; CD25
Gene ID	3559
SwissProt ID	P01589.

Application

Dilution Ratio	IHC: 1:50-1:100
Molecular Weight	-

Product Name: IL-2 Receptor alpha Mouse Monoclonal Antibody
Catalog #: AMM00727




Background

Receptor for interleukin-2.


Research Area

Immunology

Image Data



Immunohistochemistry analysis of paraffin-embedded Human liver tissue using IL-2 Receptor alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IL-2 Receptor alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.