

Product Name: Rabbit IgG Light Chain (7F9) Mouse Monoclonal Antibody
Catalog #: AMM00876

Summary

Production Name	Rabbit IgG Light Chain (7F9) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Rabbit

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	Rabbit IgG
Alternative Names	Antibody; Ig gamma chain C region; Ig kappa-b4 chain C region; Ig kappa-b5 chain C region; Ig kappa-b9 chain C region; Ig lambda chain C region; K-BAS.
Gene ID	-
SwissProt ID	-.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

**Product Name: Rabbit IgG Light Chain (7F9) Mouse
Monoclonal Antibody
Catalog #: AMM00876**




Background

Constant region of immunoglobulin light chains. Immunoglobulins, also known as antibodies, are membrane-bound or secreted glycoproteins produced by B lymphocytes.

Research Area

Immunology

Image Data



Western blot analysis of Reduced Rabbit IgG, H+L using Rabbit IgG Light Chain (7F9) antibody.

Note

For research use only.