

Product Name: Prealbumin (1F6) Mouse Monoclonal Antibody**Catalog #: AMM00893**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Antigen Information

Gene Name	TTR
Alternative Names	Prealbumin; TBPA; Amyloidosis I; Transthyretin
Gene ID	7276
SwissProt ID	P02766
Immunogen	A synthetic peptide of human Prealbumin


Background

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.


Research Area

Neuroscience


Image Data



Immunohistochemistry analysis of paraffin-embedded Human Liver Carcinoma Tissue using Prealbumin (1F6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Prealbumin (1F6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Prealbumin (1F6) in Hela, Jurkat, mouse brain, MOLT4, Human serum lysates using Prealbumin (1F6) antibody