

Product Name: PRMT6 Mouse Monoclonal Antibody
Catalog #: AMM00960



Summary

Production Name	PRMT6 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC-P,ICC/IF
Reactivity	Human,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	PRMT6
Alternative Names	HRMT1L6; PRMT 6
Gene ID	55170
SwissProt ID	Q96LA8.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa

Background

Arginine methyltransferase that can catalyze the formation of both omega-N monomethylarginine (MMA) and

Product Name: PRMT6 Mouse Monoclonal Antibody
Catalog #: AMM00960




asymmetrical dimethylarginine (aDMA), with a strong preference for the formation of aDMA. Preferentially methylates arginyl residues present in a glycine and arginine-rich domain and displays preference for monomethylated substrates.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of PRMT6 in COS7, Jurkat, C6, 3T3 and HeLa lysates using PRMT6 antibody.

Note

For research use only.