

**Product Name: Pyruvate Dehydrogenase E2 Mouse
Monoclonal Antibody
Catalog #: AMM00969**



Summary

Production Name	Pyruvate Dehydrogenase E2 Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB, ICC/IF, IP
Reactivity	Human, Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	DLAT
Alternative Names	DLAT; DLTA; E2; PBC; PDCE2
Gene ID	1737
SwissProt ID	P10515.

Application

Dilution Ratio	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
Molecular Weight	Calculated MW: 69 kDa; Observed MW: 69 kDa

**Product Name: Pyruvate Dehydrogenase E2 Mouse
Monoclonal Antibody
Catalog #: AMM00969**




Background

The pyruvate dehydrogenase complex catalyzes the overall conversion of pyruvate to acetyl-CoA and CO₂. It contains multiple copies of three enzymatic components: pyruvate dehydrogenase (E1), dihydrolipoamide acetyltransferase (E2) and lipoamide dehydrogenase (E3).

Research Area

Signal Transduction

Image Data



Western blot analysis of Pyruvate Dehydrogenase E2 in COS7, Jurkat and HeLa lysates using Pyruvate Dehydrogenase E2 antibody.

Note

For research use only.