

Product Name: Nrf1 Mouse Monoclonal Antibody**Catalog #: AMM00972**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 72 kDa

Antigen Information

Gene Name	NRF1
Alternative Names	alpha pal; alpha palindromic binding protein; Alpha palindromic-binding protein; Alpha-pal; locus control region factor 1; NFE2 related factor 1; NRF-1; Nrf1; NRF1_HUMAN; Nuclear respiratory factor 1; transcription factor 11.
Gene ID	4899
SwissProt ID	Q16656
Immunogen	A synthetic peptide of human NRF1


Background

Transcription factor that activates the expression of the EIF2S1 (EIF2-alpha) gene. Links the transcriptional modulation of key metabolic genes to cellular growth and development. Implicated in the control of nuclear genes required for respiration, heme biosynthesis, and mitochondrial DNA transcription and replication.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TCF11/NRF1 in CEM and K562 lysates using TCF11/NRF1 antibody.