

Product Name: TBL1X (4H2) Mouse Monoclonal Antibody
Catalog #: AMM03414

Summary

Production Name	TBL1X (4H2) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	TBL1X
Alternative Names	TBL1X; TBL1; F-box-like/WD repeat-containing protein TBL1X; SMAP55; Transducin beta-like protein 1X; Transducin-beta-like protein 1; X-linked
Gene ID	6907
SwissProt ID	O60907.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 58 kDa

Product Name: TBL1X (4H2) Mouse Monoclonal Antibody
Catalog #: AMM03414


Background

F-box-like protein involved in the recruitment of the ubiquitin/19S proteasome complex to nuclear receptor-regulated transcription units.


Research Area

Neuroscience

Image Data



Western blot analysis of TBL1X in SW480 lysates using TBL1X antibody.



Immunocytochemistry analysis of TBL1X (4H2) in HeLa using TBL1X antibody.

Note

For research use only.