

Product Name: HP1 alpha (3G2) Mouse Monoclonal Antibody**Catalog #: AMM03419**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,ICC/IF
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 22 kDa; Observed MW: 26 kDa

Antigen Information

Gene Name	CBX5
Alternative Names	CBX5; HP1A; Chromobox protein homolog 5; Antigen p25; Heterochromatin protein 1 homolog alpha; HP1 alpha
Gene ID	23468
SwissProt ID	P45973
Immunogen	A synthetic peptide of human HP1 alpha

Background


Component of heterochromatin that recognizes and binds histone H3 tails methylated at 'Lys-9' (H3K9me), leading to

epigenetic repression. In contrast, it is excluded from chromatin when 'Tyr-41' of histone H3 is phosphorylated (H3Y41ph). Can interact with lamin-B receptor (LBR). This interaction can contribute to the association of the heterochromatin with the inner nuclear membrane. Involved in the formation of functional kinetochore through interaction with MIS12 complex proteins.


Research Area

Epigenetics and Nuclear Signaling


Image Data



Western blot analysis of HP1 alpha in MCF-7, 293T, 3T3, C6, COS7 and mouse brain lysates using HP1 alpha antibody.



Immunocytochemistry analysis of HP1 alpha (3G2) in HeLa using HP1 alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma using HP1 alpha antibody (3G2). High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.