

**Product Name: Nucleophosmin (2D9) Mouse
Monoclonal Antibody
Catalog #: AMM03420**

Summary

Production Name	Nucleophosmin (2D9) Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ChIP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	NPM1
Alternative Names	B23; NPM
Gene ID	4869
SwissProt ID	P06748.

Application

Dilution Ratio	WB 1:500-1:1000,ChIP 1:20
Molecular Weight	Calculated MW: 33 kDa; Observed MW: 38 kDa

**Product Name: Nucleophosmin (2D9) Mouse
Monoclonal Antibody
Catalog #: AMM03420**




Background

Involved in diverse cellular processes such as ribosome biogenesis, centrosome duplication, protein chaperoning, histone assembly, cell proliferation, and regulation of tumor suppressors p53/TP53 and ARF. Binds ribosome presumably to drive ribosome nuclear export. Associated with nucleolar ribonucleoprotein structures and bind single-stranded nucleic acids. Acts as a chaperonin for the core histones H3, H2B and H4. Stimulates APEX1 endonuclease activity on apurinic/apyrimidinic (AP) double-stranded DNA but inhibits APEX1 endonuclease activity on AP single-stranded RNA. May exert a control of APEX1 endonuclease activity within nucleoli devoted to repair AP on rDNA and the removal of oxidized rRNA molecules. In concert with BRCA2, regulates centrosome duplication. Regulates centriole duplication: phosphorylation by PLK2 is able to trigger centriole replication. Negatively regulates the activation of EIF2AK2/PKR and suppresses apoptosis through inhibition of EIF2AK2/PKR autophosphorylation. Antagonizes the inhibitory effect of ATF5 on cell proliferation and relieves ATF5-induced G2/M blockade (PubMed:22528486). In complex with MYC enhances the transcription of MYC target genes (PubMed:25956029).

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of NPM1 in HeLa, Jurkat and 3T3 lysates using NPM1 antibody.

Note

For research use only.