

Product Name: GATA6 (3H3) Mouse Monoclonal Antibody
Catalog #: AMM03431

Summary

Production Name	GATA6 (3H3) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	GATA6
Alternative Names	GATA6; Transcription factor GATA-6; GATA-binding factor 6
Gene ID	2627
SwissProt ID	Q92908.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 60 kDa; Observed MW: 60 kDa

Background


Product Name: GATA6 (3H3) Mouse Monoclonal Antibody
Catalog #: AMM03431

Thought to be important for regulating terminal differentiation and/or proliferation.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of GATA6 in HeLa and 293 lysates using GATA6 antibody.

Note

For research use only.