
Product Name: Smad5 (4B10) Mouse Monoclonal Antibody**Catalog #: AMM03437**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ICC/IF,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC/IF 1:50-1:200,FC 1:50-1:100
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 60 kDa

Antigen Information

Gene Name	SMAD5 DKFZp781C1895; DKFZp781O1323; Dwfc; hSmad 5; hSmad5; JV5 1; JV5-1; MAD homolog 5; MAD mothers against decapentaplegic homolog 5; MAD; mothers against decapentaplegic
Alternative Names	homolog 5; MADH 5; MADH5; Mothers against decapentaplegic homolog 5; Mothers against DPP homolog 5; MusMLP; SMA and MAD related protein 5; SMAD 5; SMAD family member 5; SMAD mothers against DPP homolog 5; Smad5; Smad5; SMAD5_HUMAN.
Gene ID	4090
SwissProt ID	Q99717
Immunogen	A synthetic peptide corresponding to target protein


Background

Transcriptional modulator activated by BMP (bone morphogenetic proteins) type 1 receptor kinase. Smad5 is a receptor-regulated Smad (R-Smad). Smad5 is required for normal development of the cardiovascular system in vivo; lack of the Smad5 gene results in apoptosis of cardiac myocytes. 3 Upregulation of Smad5 has been reported to mediate apoptosis of gastric epithelial cells induced by Helicobacter pylori infection. Tissue specificity: Ubiquitous.


Research Area

Signal Transduction


Image Data



Western blot analysis of SMAD5 (Cterminus) in HeLa, Jurkat and K562 lysates using SMAD5 (Cterminus) antibody.



Immunocytochemistry analysis of Smad5 in HeLa cells using Smad5 antibody.



Flow Cytometry analysis of Jurkat cells stained with SMAD5 (red). Black line histogram represents the isotype control, normal mouse IgG.