

Product Name: beta Arrestin 1 (6A1) Mouse Monoclonal Antibody
Catalog #: AMM03464

Summary

Production Name	beta Arrestin 1 (6A1) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	ARRB1
Alternative Names	ARB1; ARR1; ARRB1; ARRB1_HUMAN; Arrestin 2; Arrestin beta 1; Arrestin beta-1; Beta-arrestin-1.
Gene ID	408
SwissProt ID	P49407.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 51 kDa

Product Name: beta Arrestin 1 (6A1) Mouse Monoclonal Antibody
Catalog #: AMM03464


Background

Members of arrestin/beta-arrestin protein family are thought to participate in agonist-mediated desensitization of G-protein-coupled receptors and cause specific dampening of cellular responses to stimuli such as hormones, neurotransmitters, or sensory signals. Arrestin beta 1 is a cytosolic protein and acts as a cofactor in the beta-adrenergic receptor kinase (BARK) mediated desensitization of beta-adrenergic receptors. Besides the central nervous system, it is expressed at high levels in peripheral blood leukocytes, and thus the BARK/beta-arrestin system is believed to play a major role in regulating receptor-mediated immune functions. Alternatively spliced transcripts encoding different isoforms of arrestin beta 1 have been described. [provided by RefSeq, Jan 2011]

Research Area

Signal Transduction

Image Data



Western blot analysis of betaArrestin 1 in K562 and Hela lysates using betaArrestin 1 antibody.

Note

For research use only.