

Product Name: FYB (2H9) Mouse Monoclonal Antibody
Catalog #: AMM03470



Summary

Production Name	FYB (2H9) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	FYB1
Alternative Names	ADAP; Fyb; FYN binding protein; p120/p130; SLAP130
Gene ID	2533
SwissProt ID	O15117.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 85 kDa; Observed MW: 120 kDa

Background

Acts as an adapter protein of the FYN and LCP2 signaling cascades in T cells. Modulates the expression of interleukin-2 (IL-

Product Name: FYB (2H9) Mouse Monoclonal Antibody
Catalog #: AMM03470




2). Involved in platelet activation. Prevents the degradation of SKAP1 and SKAP2.

Research Area

Signal Transduction

Image Data



Western blot analysis of FYB in 10, 15 and 20ug Jurkat lysates using FYB antibody.

Note

For research use only.