

Product Name: Phostensin (4B12) Mouse Monoclonal Antibody
Catalog #: AMM03494



Summary

Production Name	Phostensin (4B12) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	PPP1R18
Alternative Names	HKMT1098; KIAA1949.
Gene ID	170954
SwissProt ID	Q6NYC8.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 68 kDa; Observed MW: 110, 26 kDa

Background

Product Name: Phostensin (4B12) Mouse Monoclonal Antibody
Catalog #: AMM03494




May target protein phosphatase 1 to F-actin cytoskeleton.

Research Area

Signal Transduction

Image Data



Western blot analysis of Phostensin in Hela, Raji, Ramos and Leucocyte lysates using Phostensin antibody.

Note

For research use only.