

Product Name: SMYD4 (7C8) Mouse Monoclonal Antibody**Catalog #: AMM03508**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 89 kDa; Observed MW: 89 kDa

Antigen Information

Gene Name	SMYD4
Alternative Names	ZMYND21.
Gene ID	114826
SwissProt ID	Q8IYR2
Immunogen	


Background

Plays a critical role in cardiac development.


Research Area

Cancer

Image Data



Western blot analysis of SMYD4 in Hela, MCF-7 and A431 lysates using SMYD4 antibody.



Western blot analysis of SMYD4 (7C8) in Hela, MCF-7 and MDA-MB-468 lysates using SMYD4 antibody.