

Product Name: BRAP (1E7) Mouse Monoclonal Antibody
Catalog #: AMM03514

Summary

Production Name	BRAP (1E7) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	BRAP
Alternative Names	IMP; BRAP2; RNF52.
Gene ID	8315
SwissProt ID	Q7Z569.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 67 kDa; Observed MW: 67 kDa

Background

Product Name: BRAP (1E7) Mouse Monoclonal Antibody
Catalog #: AMM03514




Negatively regulates MAP kinase activation by limiting the formation of Raf/MEK complexes probably by inactivation of the KSR1 scaffold protein.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of BRAP (1E7) in CHO-K1 lysates, transfected with a Human pFLAGCMV2BRAP construct (A) or transfected with a Human pFLAGCMV2 construct (B), using BRAP antibody.



Western blot analysis of BRAP (1E7) in A549 and HUVEC lysates using BRAP antibody.

Note

For research use only.