

Product Name: GRAP2 (5C7) Mouse Monoclonal Antibody
Catalog #: AMM03541

Summary

Production Name	GRAP2 (5C7) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	GRAP2
Alternative Names	P38; GADS; GRID; GRPL; GrbX; Mona; GRB2L; GRBLG; Grf40; GRAP-2
Gene ID	9402
SwissProt ID	O75791.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Background


Product Name: GRAP2 (5C7) Mouse Monoclonal Antibody
Catalog #: AMM03541

Interacts with SLP-76 to regulate NF-AT activation. Binds to tyrosine-phosphorylated shc.

Research Area

Signal Transduction

Image Data



Western blot analysis of GRAP2 in GRAP2 positive (CEM, Molt4, Jurkat) and negative (3T3, K562, Hela) lysates using GRAP2 antibody.

Note

For research use only.