

Product Name: Lamin B2 (8F6) Mouse Monoclonal Antibody
Catalog #: AMM03554

Summary

Production Name	Lamin B2 (8F6) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	LMNB2
Alternative Names	LMN2; LAMB2
Gene ID	84823
SwissProt ID	Q03252.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 70 kDa; Observed MW: 70 kDa

Product Name: Lamin B2 (8F6) Mouse Monoclonal Antibody
Catalog #: AMM03554




Background

Lamins are components of the nuclear lamina, a fibrous layer on the nucleoplasmic side of the inner nuclear membrane, which is thought to provide a framework for the nuclear envelope and may also interact with chromatin.


Research Area

Cell Biology

Image Data



Western blot analysis of Lamin B2 in HCT116, mouse Embryo and Jurkat lysates using Lamin B2 antibody.



Immunoprecipitation analysis of Lamin B2 (8F6) in HeLa lysates using Lamin B2 antibody.

Note

For research use only.