
Product Name: SSB (4B10) Mouse Monoclonal Antibody**Catalog #: AMM03555**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 50 kDa

Antigen Information

Gene Name	SSB
Alternative Names	La; LARP3; La/SSB
Gene ID	6741
SwissProt ID	P05455
Immunogen	


Background

Binds to the 3' poly(U) termini of nascent RNA polymerase III transcripts, protecting them from exonuclease digestion and facilitating their folding and maturation.

Research Area

Immunology

Image Data



Western blot analysis of La Antigen in MCF-7, A549, Jurkat and HeLa lysates using La Antigen antibody.