

Product Name: HARS (4H5) Mouse Monoclonal Antibody
Catalog #: AMM03557

Summary

Production Name	HARS (4H5) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human,Mouse,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	HARS1
Alternative Names	HRS; CMT2W; USH3B
Gene ID	3035
SwissProt ID	P12081.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 57 kDa; Observed MW: 57 kDa

Background


Product Name: HARS (4H5) Mouse Monoclonal Antibody
Catalog #: AMM03557

Cytoplasmic histidine--tRNA ligase (Probable). Plays a role in axon guidance.

Research Area

Immunology

Image Data



Western blot analysis of HARS in T47D, 3T3, COS7, C6 and HeLa lysates using HARS antibody.

Note

For research use only.