
Product Name: EpCAM (1D5) Mouse Monoclonal Antibody**Catalog #: AMM03641**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 39 kDa

Antigen Information

Gene Name	EPCAM EPCAM; GA733-2; M1S2; M4S1; MIC18; TACSTD1; TROP1; Epithelial cell adhesion molecule; Ep-CAM; Adenocarcinoma-associated antigen; Cell surface glycoprotein Trop-1; Epithelial
Alternative Names	cell surface antigen; Epithelial glycoprotein; EGP; Epithelial glycoprotein 314; EGP314; hEGP314; KS 1/4 antigen; KSA; Major gastrointestinal tumor-associated protein GA733-2; Tumor-associated calcium signal transducer 1; CD326
Gene ID	4072
SwissProt ID	P16422
Immunogen	A synthetic peptide of human EpCAM


Background

May act as a physical homophilic interaction molecule between intestinal epithelial cells (IECs) and intraepithelial lymphocytes (IELs) at the mucosal epithelium for providing immunological barrier as a first line of defense against mucosal infection. Plays a role in embryonic stem cells proliferation and differentiation. Up-regulates the expression of FABP5, MYC and cyclins A and E.


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of EpCAM (1D5) in CHO-K1 transfected by FlagEpCAM, MCF-7 (EpCAM positive), HCT116 (EpCAM positive), and HeLa (EpCAM negative) lysates using EpCAM (1D5) antibody.



Immunoprecipitation analysis of EpCAM (1D5) in HCT116 and HeLa lysates using EpCAM antibody.