

Product Name: ICAM1 (4F5) Mouse Monoclonal Antibody
Catalog #: AMM03648

Summary

Production Name	ICAM1 (4F5) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	ICAM1
Alternative Names	ICAM1; Intercellular adhesion molecule 1; ICAM-1; Major group rhinovirus receptor; CD antigen CD54
Gene ID	3383
SwissProt ID	P05362.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 58 kDa; Observed MW: 96 kDa

Product Name: ICAM1 (4F5) Mouse Monoclonal Antibody
Catalog #: AMM03648


Background

ICAM1 a type I membrane protein of the immunoglobulin superfamily. Is a ligand for the leukocyte adhesion LFA-1 protein (Integrin alpha-L/beta-2) and a Rhinovirus receptor. Typically expressed on endothelial cells and cells of the immune system.


Research Area

Immunology

Image Data



Western blot analysis of CD54(ICAM1) in SW480 and Raji lysates using CD54(ICAM1) antibody.



Immunocytochemistry analysis of ICAM1 (4F5) in SW480 using CD54(ICAM1) antibody.

Note

For research use only.