

Product Name: Hexokinase II (3B9) Mouse Monoclonal Antibody
Catalog #: AMM03657



Summary

Production Name	Hexokinase II (3B9) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	HK2
Alternative Names	HKII; HXK2
Gene ID	3099
SwissProt ID	P52789.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 102 kDa; Observed MW: 102 kDa

Background

Product Name: Hexokinase II (3B9) Mouse Monoclonal Antibody
Catalog #: AMM03657




Miscellaneous In vertebrates there are four major glucose-phosphorylating isoenzymes, designated hexokinase I, II, III and IV (glucokinase).

Research Area

Signal Transduction

Image Data



Western blot analysis of Hexokinase II in COS7, K562 and HeLa lysates using Hexokinase II antibody.

Note

For research use only.