

Product Name: Lyn (1C8) Mouse Monoclonal Antibody
Catalog #: AMM03690



Summary

Production Name	Lyn (1C8) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	LYN
Alternative Names	LYN; JTK8; Tyrosine-protein kinase Lyn; Lck/Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn
Gene ID	4067
SwissProt ID	P07948.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 59 kDa; Observed MW: 56,53 kDa

Background

Product Name: Lyn (1C8) Mouse Monoclonal Antibody
Catalog #: AMM03690




Lyn a tyrosine kinase of the Src family that plays a role in B cell antigen-receptor signaling.

Research Area

Signal Transduction

Image Data



Western blot analysis of Lyn in Ramos, Raji, K562 and Hela lysates using Lyn antibody.

Note

For research use only.