

Product Name: OVA Mouse Monoclonal Antibody
Catalog #: AMM80563



Summary

Production Name	OVA Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,ELISA
Reactivity	Other

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	OVA
Alternative Names	OVA
Gene ID	
SwissProt ID	. Ovalbumin

Application

Dilution Ratio	WB:1:500-1:2000,ELISA:1:10000
Molecular Weight	45kDa

Background

Ovalbumin (OVA) is used in variety of laboratories products which contain OVA as a carrier protein. To raise antibodies against small non-immunogenic molecules, called haptens, it is often necessary to couple them to carrier protein, such as


Product Name: OVA Mouse Monoclonal Antibody
Catalog #: AMM80563



OVA. This antibody is a good reagent that may be used for the development of detection assay in ELISA or western blot to test the remaining OVA for clinical products in Vivo or research products In Vitro.

Research Area

Image Data



Western blot analysis using OVA mouse mAb against OVA protein.

Note

For research use only.