

Product Name: FRK Mouse Monoclonal Antibody
Catalog #: AMM80667



Summary

Production Name	FRK Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	FRK
Alternative Names	GTK; RAK; PTK5
Gene ID	2444.0
SwissProt ID	P42685 . Purified recombinant fragment of FRK (aa2-300) expressed in E. Coli.

Application

Dilution Ratio	IHC:1:200-1:1000,ELISA:1:10000
Molecular Weight	/

Background


FRK: fyn-related kinase. The protein encoded by this gene belongs to the TYR family of protein kinases. This tyrosine kinase is a nuclear protein and may function during G1 and S phase of the cell cycle and suppress growth.

Product Name: FRK Mouse Monoclonal Antibody
Catalog #: AMM80667



Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human Breast tissues using FRK mouse mAb.

Note

For research use only.