

Product Name: PAR1 Mouse Monoclonal Antibody**Catalog #: AMM80681**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	/

Antigen Information


Gene Name	PAR1
Alternative Names	PAR1; D15S227E
Gene ID	145624.0
SwissProt ID	
Immunogen	Purified recombinant fragment of PAR1 expressed in E. Coli.

Background


PAR1: Prader-Willi/Angelman region-1 .

Research Area

Image Data



Flow cytometric analysis of hela cells using PAR1 mouse mAb (green) and negative control (red).



Flow cytometric analysis of k562 cells using PAR1 mouse mAb (green) and negative control (red).