
Product Name: human IgG (Fc Specific) Mouse Monoclonal Antibody**Catalog #: AMM80731**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	/

Antigen Information


Gene Name	human IgG (Fc Specific)
Alternative Names	human IgG (Fc Specific)
Gene ID	
SwissProt ID	/
Immunogen	Purified fusion protein with human IgG (Fc Specific) tag.

Background

Specificity of antibody has been tested using Western blot, immunofluorescence, and ELISA. The product reacts with human IgG (Fc) with no significant reactivity with human IgG (Fab2), IgM or other serum proteins. We didn't test crossing with other species IgG (Fc).

Research Area

Image Data



Western blot analysis using human IgG (Fc specific) mouse mAb against human serum (1).