
Product Name: NTRK3 Mouse Monoclonal Antibody**Catalog #: AMM80770**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000
Molecular Weight	94kDa

Antigen Information

Gene Name	NTRK3
Alternative Names	TRKC; gp145(trkC)
Gene ID	4916.0
SwissProt ID	Q16288
Immunogen	Purified recombinant extracellular fragment of human NTRK3 (aa32-429) fused with hlgGfc tag expressed in HEK293 cells.


Background

NTRK3 (neurotrophic tyrosine kinase, receptor, type 3), it is a member of the neurotrophic tyrosine receptor kinase (NTRK) family and plays an important role in the development and maintenance of neural tissues. This kinase is a membrane-bound

receptor that, upon neurotrophin binding, phosphorylates itself and members of the MAPK pathway. Signalling through this kinase leads to cell differentiation and may play a role in the development of proprioceptive neurons that sense body position. Mutations in this gene have been associated with medulloblastomas, secretory breast carcinomas and other cancers.

Research Area

Image Data



Western blot analysis using NTRK3 mouse mAb against extracellular domain of human NTRK3 (aa32-429).