
Product Name: NGFR Mouse Monoclonal Antibody**Catalog #: AMM80792**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ICC,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	45kDa

Antigen Information

Gene Name	NGFR
Alternative Names	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR; NGFR
Gene ID	4804.0
SwissProt ID	P08138
Immunogen	Purified recombinant fragment of human NGFR expressed in E. Coli.

Background


Nerve growth factor receptor contains an extracellular domain containing four 40-amino acid repeats with 6 cysteine residues at conserved positions followed by a serine/threonine-rich region, a single transmembrane domain, and a 155-amino acid cytoplasmic domain. The cysteine-rich region contains the nerve growth factor binding domain. NGFR p75 plays a central role

in the regulation of cell number by apoptosis in the developing CNS. During early development, activation of NGFR p75 by NGF induces apoptotic cell death in some neuronal cells, probably through activation of the sphingomyelinase/ceramide pathway, the ICE like proteases and the JNK pathway. In rat Schwann cells, NGF binding to NGFR p75 activates NF kappaB, possibly to modulate Schwann cell migration during nerve regeneration. CD271 has recently been described as being expressed in mesenchymal stem cells (bone marrow stromal cells).


Research Area

Apoptosis, TGF-beta signaling pathway, PI3K-Akt signaling pathway, Hippo signaling pathway

Image Data



Immunofluorescence analysis of EC cells using NGFR mouse mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of EC cells using NGFR mouse mAb (green) and negative control (purple).