
Product Name: KLHL21 Mouse Monoclonal Antibody**Catalog #: AMM80857**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:200-1:1000,ELISA 1:5000-1:20000
Molecular Weight	66kDa

Antigen Information


Gene Name	KLHL21
Alternative Names	KIAA0469; MGC99635; KLHL21
Gene ID	9903.0
SwissProt ID	Q9UJP4
Immunogen	Purified recombinant fragment of human KLHL21 expressed in E. Coli.

Background

Probable substrate-specific adapter of an E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human brain tissues using KLHL21 mouse mAb with DAB staining.