

Product Name: KLHL11 Mouse Monoclonal Antibody
Catalog #: AMM80862



Summary

Production Name	KLHL11 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,ICC,FC,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	KLHL11
Alternative Names	FLJ10572; KLHL11
Gene ID	55175.0
SwissProt ID	Q9NVR0 .Purified recombinant fragment of human KLHL11 expressed in E. Coli.

Application

Dilution Ratio	WB:1:500-1:2000,ICC:1:200-1:1000,FC:1:200-1:400,ELISA:1:10000
Molecular Weight	80kDa

Background


Kelch-like protein 11, Subcellular location Secreted.

Product Name: KLHL11 Mouse Monoclonal Antibody
Catalog #: AMM80862




Research Area


Image Data



Western blot analysis using KLHL11 mouse mAb against Hela (1) and MCF-7 (2) cell lysate.



Immunofluorescence analysis of 3T3-L1 cells using KLHL11 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of K562 cells using KLHL11 mouse mAb (green) and negative control (purple).

Product Name: KLHL11 Mouse Monoclonal Antibody
Catalog #: AMM80862



Note

For research use only.