
Product Name: CHUK Mouse Monoclonal Antibody**Catalog #: AMM80891**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	85kDa

Antigen Information

Gene Name	CHUK
Alternative Names	IKK1; IKKA; IKBKA; TCF16; NFKB1A; IKK-alpha; CHUK
Gene ID	1147.0
SwissProt ID	O15111
Immunogen	Purified recombinant fragment of human CHUK expressed in E. Coli.


Background

This gene encodes a member of the serine/threonine protein kinase family. The encoded protein, a component of a cytokine-activated protein complex that is an inhibitor of the essential transcription factor NF-kappa-B complex, phosphorylates sites that trigger the degradation of the inhibitor via the ubiquitination pathway, thereby activating the transcription factor.

Research Area

PI3K-Akt signaling pathway, MAPK signaling pathway

Image Data



Western blot analysis using CHUK mouse mAb against Raji (1), Jurkat (2) and THP-1 (3) cell lysate.