

---

**Product Name: GKAP Mouse Monoclonal Antibody****Catalog #: AMM81013**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	IHC,ICC,ELISA
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	PBS containing 0.03% sodium azide.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000
<b>Molecular Weight</b>	109kDa

**Antigen Information**


<b>Gene Name</b>	GKAP
<b>Alternative Names</b>	GKAP; DLGAP1; DAP-1; hGKAP; SAPAP1; FLJ38442; MGC88156; DAP-1-BETA; DAP-1-ALPHA
<b>Gene ID</b>	9229.0
<b>SwissProt ID</b>	O14490
<b>Immunogen</b>	Purified recombinant fragment of human GKAP expressed in E. Coli.

**Background**


Function: Part of the postsynaptic scaffold in neuronal cells.Tissue specificity: Expressed in brain.

**Research Area**


## Image Data




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Immunofluorescence analysis of NIH/3T3 cells using GKAP mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using GKAP mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophagus tissues using GKAP mouse mAb with DAB staining.